

Red-Light Hip Heating Pad



① Hip Heating Pad ② Product Strap ③ Charging Cable

Product Manual

Please read this manual carefully before use and keep it properly for future reference

Product Introduction

The Red-light heated hip pad is crafted from high – quality diving fabric and features an overheat protection design, ensuring a comfortable and safe user experience. It is ergonomically designed and comes with a long belt for added comfort. The heating therapy can relieve muscle pain, relax muscles, and reduce stress, especially effective in alleviating pain in the hip joints and thighs.

Notice:

The heated hip pad can help relieve pain, yet it does not have therapeutic properties

Product Parameters

Red light LED beads	Output Specifications	Battery Capacity	Heating wire	Full power range
105 pieces	5V/3A	5000mAh	Carbon fiber material	Approximately 50 min
Operating Temperature				
High Mode	Medium Mode		Low Mode	
150°F	130°F		110°F	

Charge the Product

The product is equipped with a built-in 5000mAh battery




Low Battery Reminder



Full Charge Reminder

How to Use



1. Please fully charge the product prior to use. Battery life is approximately 50 minutes at the highest power setting.
2. Place the Hip Heating Pad on your left or right hip, then adjust the straps for a comfortable, secure fit.
3. Adjust the waist and leg Velcro straps to ensure the vibration motors fit snugly against your skin.
4. Press and hold the Temperature & Vibration button to activate the device.
5. Press the button briefly to cycle through the 3 adjustable modes. The device will automatically shut off after 25 minutes of continuous use in any mode. To continue using it, simply turn it on again.

Cleaning and Maintenance

Before cleaning, power off the device by pressing and holding the function button until the indicator light goes out. Gently wipe the heating pad with a neutral cleaner (keep the controller dry), then lay flat to air-dry. Avoid harsh detergents and stiff brushes.

1. Do not clean the Hip Heating Pad in a washing machine.
2. Do not dry clean, as this may damage the heating element. Wiping with water below 30°C (86°F) is recommended.
3. Do not immerse the controller in water. If the controller surface becomes dirty, wipe it with a soft cloth and neutral cleaner.
4. Make sure the device is completely cooled before storing.
Store the product flat in a clean, dry place.



Not Washable in Washing Machine



Do not wring out



Do not tumble dry



Do not iron.

Attention

1. When using, wear thin clothing such as yoga wear. Do not wear thick clothing, as it will reduce heating efficiency.
2. Do not wear this product directly against bare skin (without clothing as a barrier), as it may cause burns.
3. Before first use, cycle through high, medium, and low temperatures to test your skin's heat sensitivity, then select a suitable temperature.
4. After use, fully power off the device by pressing and holding the function button until the indicator light turns off.
5. When storing, keep the product in a cool, dry environment, away from high temperatures, direct sunlight (e.g., inside a car, near radiators), and humid areas.
6. Do not fold or compress the heating pad forcefully. Store it in a well-ventilated, flat area and keep it away from flammable materials.
7. If not in use for an extended period, maintain the battery charge at approximately 50% and recharge it every 3 months to prevent battery depletion and failure.
8. This product is prohibited for use by pregnant women, young children, and individuals who cannot communicate their needs.